

Wound dressing for moist wound healing


Plus moist™

Plus moist is a new type of wound dressing developed on moist wound healing concept. Plus moist gives superior non-adherent ability in a moist environment, and does not form a scab. Because Plus moist uses self-regulating absorption of the exudate fluid, Plus moist provides moist environment without rashes or maceration (skin deterioration) to various wounds. Plus moist is ideal for traumatic wounds (Ex. burns, abrasions, crushed wounds, suture wounds etc.), chronic wounds (Ex. ulcers, pressure ulcers) and dermatitis (Ex. rashes, eczema, contagious impetigo).

Features:

- 1) Non-adherent (It does not dry wound surface).
- 2) Superior absorption ability.
- 3) Unique Self-Regulating Absorption of wound exudate.
- 4) No rashes or maceration (skin deterioration).
- 5) Able to cut for any size and shape.

Direction for use:



White mesh film side of dressing to wound.
Able to cut for any size and shape.





Fig. Function as a self-regulating absorption of the exudate



Main Product Line

Code	Product Name	Size	Package	Remark
V1	Plus moist V	250 × 200mm	10pcs / bag	
P1	Plus moist P	250 × 200mm	3pcs / bag	
P3		125 × 125mm		
W1	Plus moist W	250 × 200mm	10pce / box	EB sterilization
W2		125 × 200mm		
W3		125 × 125mm		

WARNING

Do not reuse. Single use only. Do not use on infected wound which indicates signs such as fever, pain, redness, itching, rash or burning etc. Diabetics and those with poor circulation or wound healing problems should not use this product.



【Manufacturer】
ZUIKO MEDICAL CORPORATION
 4-3-50, Torikaikami, Settsu, Osaka, Japan
 Tel: 81-72-653-8877 Fax: 81-72-653-8876
 Mail: plus_moist@zuiko.co.jp
 Website: www.zuiko-medical.co.jp/